Stabilized HEB-QCL heterodyne spectrometer at superterahertz

Y. Ren^{1,2,3,*,} J.N. Hovenier¹, D.J. Hayton⁴, M. Cui⁴, J.R. Gao^{1,4*,} T.M. Klapwijk¹, S.C. Shi², T-Y. Kao⁵, Q. Hu⁵, and J. L. Reno⁶

1Kavli Institute of NanoScience, Delft University of Technology, Lorentzweg 1, 2628 CJ Delft, The Netherlands 2Purple Mountain Observatory (PMO), Chinese Academy of Sciences, 2 West Beijing Road, Nanjing, JiangSu 210008, China.

3Graduate School, Chinese Academy of Sciences, 19A Yu Quan Road, Beijing 100049, China 4SRON Netherlands Institute for Space Research, Sorbonnelaan 2, 3584 CA Utrecht, The Netherlands 5Department of Electrical Engineering and Computer Science, Massachusetts Institute of Technology, Cambridge, Massachusetts 02139, USA. 6Sandia National Laboratories, Albuquerque, NM 87185-0601, U.S.A.

bSandia National Laboratories, Albuquerque, NM 8/185-0601, U.S.A * Contact: y.ren@tudelft.nl

Abstract—In the terahertz (THz) frequency range a high resolution heterodyne spectrometer is of crucial importance for astronomical observations and atmospheric remote sensing applications, based on its combination of high spectral resolution and sensitivity. As the mixer, a superconducting NbN hot-electron bolometer (HEB) mixer has demonstrated excellent sensitivities up to 5.3 THz. As the local oscillator (LO), terahertz quantum cascade lasers (QCLs) have shown advantages at frequencies above 2 THz, based on their single mode emission, wide frequency operating range, high output power and long term stability. Several progresses have been made for a THz QCL to be used as the local oscillator (LO) in a heterodyne receiver, such as a heterodyne spectroscopy measurement in the lab.[1]

Here we report a new experiment on high-resolution heterodyne spectrometer using a 3.5 THz QCL as LO and a HEB as mixer by stabilizing both frequency and amplitude of the QCL. We have already achieved the following results. The frequency locking of the QCL was demonstrated by using a molecular absorption line, a proportional-integral-derivative (PID) controller, and a direct power detector [2]. The intensity of the QCL emission is also stabilized by means of swing-arm actuator placed in the beam path using a second PID controller [3].

Y. Ren, J.N. Hovenier, R. Higgins, J.R. Gao, T.M. Klapwijk, S.C. Shi, B. Klein, T-Y. Kao, Q. Hu, and J. L. Reno, High-resolution heterodyne spectroscopy using a tunable quantum cascade laser around 3.5 THz. Appl. Phys. Lett., **98**, 231109(2011).
Y. Ren, J.N. Hovenier, M. Cui, D.J. Hayton, J.R. Gao, T.M. Klapwijk, S.C. Shi, T-Y. Kao, Q. Hu, and J. L. Reno, Frequency locking of 3.5 THz quantum cascade lasers using a gas cell, Appl. Phys. Lett. (submitted)
D.J. Hayton, J.R. Gao, J.W. Kooi, Y. Ren, W. Zhang and G. de Lange, Stabilized hot electron bolometer heterodyne receiver

[3] D.J. Hayton, J.R. Gao, J.W. Kooi, Y. Ren, W. Zhang and G. de Lange, Stabilized hot electron bolometer heterodyne receiver at 2.5 THz, Appl. Phys. Lett. (submitted)